

# QUICK START GUIDE



## POWER UP

Turn on the following:

- the computer and the system monitor
- once computer desktop screen is on, turn the **KEY** to the **ON** position on the front of the CoreVision

## START SOFTWARE

Double click on the **BIOPTICS** icon in the centre of the monitor

## CALIBRATE

Calibration will occur automatically once you enter the software - takes about 3 minutes

## PATIENT DATA

- Select:**
- Start Procedure
  - Patient Data
  - Enter manually or select from Worklist
  - Your Worklist Computer Icon
  - Connect
  - Find Worklist
  - OK - Search MWL
  - Locate patient's name with proper exam and double click over it
  - OK

## X-RAY IMAGING

- Place core specimens on tray provided
- Open, and close, tray door by pressing the **black button**
- Press and release the **green button** to x-ray cores - they will appear on screen after about 5 seconds

## IMAGE MANIPULATION

- Hold down the **left mouse button** and drag the mouse to change brightness and contrast for window/level function
- Click on **Enhance Image** to add more contrast to califications
- **Full Resolution** will enlarge the image for viewing at its full (native) x-ray resolution
- **Invert** will transpose the black and white colors for positive / negative film effect

## STORING TO PACS

- Choose **PACS Store** icon on toolbar and click on **PACS Storage Server** icon
- Connect
- Store
- Computer icon will stop flashing when send is complete

## PRINTING

- Choose **PACS Print** icon on toolbar and click on **PRINT Server** icon
- Connect
- Print
- Computer icon will stop flashing when print is complete

## SHUT-DOWN

- Click **EXIT**, select **SHUTDOWN** and then click OK
- Once computer monitor is blank turn the **KEY** on CoreVision to **OFF** position