

QUICK START GUIDE



POWER UP

Turn on the following switches:

- Press purple **ON/OFF** button until audible beep is heard

START SOFTWARE

Check the X-Ray chamber to ensure that it is completely empty
Double click on the **BIOPTICS** icon in the center of the monitor
Enter the password, if prompted

CALIBRATE

Calibration will occur automatically once you enter the software - takes about 3 minutes

PATIENT DATA

- Select:**
- Start Procedure
 - Patient Data
 - Enter manually or select from Worklist
 - Your Worklist Computer Icon
 - Connect
 - Find Worklist
 - OK - Search MWL
 - Locate patient's name and double click over it
 - OK

X-RAY IMAGING

- Place specimens in the X-Ray chamber, use magnification tray if needed
- Press green button to activate x-ray - image appears after a few seconds

IMAGE MANIPULATION

- Hold down the **left mouse button** and drag the mouse to change brightness and contrast for window/level function
- Click on **Enhance Image** to add more contrast to calcifications
- **Full Resolution** will enlarge the image for viewing at its full (native) x-ray resolution
- **Invert** will transpose the black and white colors for positive / negative film effect

STORING TO PACS

- Choose **PACS Store** icon on left of screen and click on **PACS Computer** icon
- Connect
- DICOM Store Selection
- Computer icon will stop flashing when send is complete

PRINTING

- Choose **PRINT** icon on left of screen and click on **PRINTER Computer** icon
- Connect
- Print DICOM Print Selection
- Computer icon will stop flashing when print is complete

SHUT-DOWN

- Click **EXIT**, select **SHUTDOWN** and then click OK
- Once computer monitor is blank press purple **ON/OFF** button